
Reveal[®]

for Multi-Treenut

INTRODUCTION AND INTENDED USE

Reveal[®] for Multi-Treenut is uniquely designed to screen environmental swabs, rinses and validated food products for the presence of any one or a combination of almond, hazelnut, pecan, walnut, cashew and/or pistachio.

Reveal for Multi-Treenut is intended for use only in an industrial food manufacturing/preparation context. Because of the problems of adequate sampling and extraction of treenuts, it is not suitable for the testing of foods to be consumed in the home or in a restaurant by allergic individuals. The test detects significant (low parts per million, or ppm) levels of almond, hazelnut, pecan, walnut, cashew and/or pistachio in clean-in-place rinses and environmental swabs.

LIMIT OF DETECTION

The test utilizes highly specific antibodies to detect almond, hazelnut, pecan, walnut, cashew and/or pistachio proteins. These treenuts can be detected from various surfaces using the provided environmental swabs or from rinses.

Utilizing the environmental swabs supplied, levels of 5 µg/100 cm² of walnut extract, 10 µg/100 cm² of pecan or hazelnut extracts and 20 µg/100 cm² of almond, cashew or pistachio extracts can be detected from Teflon surfaces. For further information on validated surfaces, please see the Reveal for Multi-Treenut validation report.

When analyzing rinses, almond, hazelnut, pecan, walnut, cashew and/or pistachio residues are detectable at a level of 5–10 ppm. The presence of cleaners and sanitizers can affect the limit of detection (LOD) in rinses.

CROSS REACTIVITY

The test was used to analyze a panel of potentially cross reacting commodities including nuts, seeds, legumes and grains. Of these, significant amounts of linseed, sunflower seed, pumpkin seed, coconut and sesame flour gave cross reactivity with the test. Approximate levels where cross reactivity was detected were 0.1% ground sunflower seed or linseed, 0.25% desiccated coconut and 1% sesame flour or ground pumpkin seed. Further information is available in the Reveal for Multi-Treenut validation report, which is available on request.

TEST PERFORMANCE

Reveal for Multi-Treenut has undergone rigorous validation to evaluate the specificity, sensitivity, robustness and intra- and inter-batch variability of the test method on rinses and environmental swabs. A validation report is available on request.

SAMPLE COMPATIBILITY

Reveal for Multi-Treenut is designed to detect low level contamination of almond, hazelnut, pecan, walnut, cashew and/or pistachio on environmental surfaces and in rinses. Although every effort has been made to validate as many variables as possible, there may be some sample types that are not suitable for testing.

The test is not suitable for testing food commodities or ingredients that contain high concentrations of the listed treenuts. At levels greater than 20% of each individual treenut in a sample the device may overload, which may result in a false negative reading on the device.

The validation of certain food ingredients may be applicable for the Reveal test format. However, it is recommended that a more broadly inclusive test method, such as Veratox or external laboratory services, be used for food ingredients. Commodity validations should be revalidated when suppliers or the manufacturing process has been changed. Validate detection by testing a positive control from the source of the allergen or ingredient that is available in your manufacturing environment and that may pose an allergen contamination risk to ensure the test can detect the allergen of concern. Please contact a Neogen representative or distributor for additional details.

MATERIALS PROVIDED

The Reveal for Multi-Treenut test pack contains the following:

1. An instruction leaflet
2. 1 foil pouch, containing 10 lilac Reveal for Multi-Treenut devices
3. 10 sachets containing Type 11 Extraction Buffer
4. 10 sample tubes and caps
5. 10 individually packaged, sterile swabs with break-off tips
6. 1 bottle of swab wetting solution

SAMPLING TECHNIQUE

1. Clean-in-place rinse water

As only a small amount of material is required for the Reveal for Multi-Treenut test, it is important to test a representative portion of liquid.

2. Environmental swabbing

The swabs supplied are intended to be used for the collection of environmental samples from which the presence of almond, hazelnut, pecan, walnut, cashew and/or pistachio can be tested. This method can be used to validate the adequacy of cleaning and/or to identify problem areas (e.g., unwanted build-up of treenuts in processing equipment).

3. Food samples

Contact a Neogen representative or distributor for additional details.

PRECAUTIONS

1. For environmental testing in an industrial food manufacturing/preparation, or labeling enforcement context only.
2. Do not use any part of the test beyond the expiry date.
3. Do not open the foil pouch until just before use.
4. Ensure the foil bag is tightly sealed after removal of a device.
5. Always store the kit between 2–8°C (35–46°F). Avoid freezing.
6. Bring kit to room temperature 18–30°C (64–86°F) prior to use.
7. The test devices use a pale green loading dye at T (test) and C (control) positions for quality purposes. The pale green lines should be visible on the devices prior to use. If the lines are not visible, please do not use the device and contact a Neogen representative for assistance.

TEST STORAGE / SAMPLE STABILITY

Store the Reveal for Multi-Treenut kit between 2–8°C (35–46°F), and use within the expiry date stated on the outer label. Extracted samples should be used within **3 hours** of extraction.

LIMITATIONS



A negative test from a surface or rinse cannot exclude the possibility that the food or swabbed surface has treenuts present since it may be distributed unevenly on the surface or in the rinse and may be below the limit of detection of the test.

Reveal for Multi-Treenut is qualitative and should only be used as a preliminary screen for the presence of almond, hazelnut, pecan, walnut, cashew and/or pistachio. The validity of results obtained with the test should preferably be viewed in conjunction with data from an alternative method such as Veratox or through external laboratory services.

RINSE SAMPLING

Remove the following and allow to equilibrate at room temperature before use (**20–30 minutes** out of refrigerator):


1. 1 sample tube and cap
 2. 1 Type 11 Extraction Buffer sachet
 3. 1 Reveal for Multi-Treenut device (in foil pouch)
1. Carefully tear/cut off the top of the Type 11 Extraction Buffer sachet and add the entire contents to the sample tube.
 2. Add the 1 mL sample to the sample tube. If a pipette is unavailable, 1 mL may be crudely estimated by half-filling a white sample tube cap. **NOTE:** Test results will be less accurate.
 3. Secure the white cap and shake for **1 minute**.

SWAB SAMPLING


Remove the following and allow to equilibrate at room temperature before use (20–30 minutes out of refrigerator):

1. 1 bottle of swab wetting solution
2. 1 sample tube and cap
3. 1 Type 11 Extraction Buffer sachet
4. 1 Reveal for Multi-Treenut device (in foil pouch)
5. 1 sterile swab


1. Carefully tear/cut off the top of the Type 11 Extraction Buffer sachet and add the entire contents to the sample tube.

2. Estimate a swabbing area of approximately 10 cm x 10 cm. Alternatively, use the swab to collect samples of contamination from problem areas (e.g., of processing equipment).


3. Gather the sample with the swab, using one of the following methods:


For dry surfaces: Remove a sterile swab from the packaging and wet with 2 drops of swab wetting solution. Swab a 10 x 10 cm area by using a crosshatch technique revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing.

For wet surfaces: Remove a sterile swab from the packaging and swab a 10 x 10 cm area by using a crosshatch technique, revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing.


4. Place the swab to the Extraction Buffer in the sample tube and carefully break off the moistened end at the prescored mark so that it remains in the tube.
5. Secure the cap of the sample tube, taking care to ensure that the stem does not prevent the tube from being properly sealed. Shake for **1 minute**.

SAMPLE TESTING

1. Remove the lid and fill it with the liquid from the tube. Any froth should remain in the tube.
2. Dip the head of the Reveal for Multi-Treenut device into the liquid in the lid. Ensure that the cavity is saturated with the liquid.
3. Leave the cavity saturated until liquid is seen running in the test window.
4. Place device on a flat surface and allow test to develop for **10 minutes**.


NOTE: It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device. Additionally, the devices are pre-stripped with pale green loading dye in positions T (test) and C (control). The loading dye assists with quality and manufacturing checks and does not impact test performance. The loading dye is removed from the test window as the sample flows through the device.

READING RESULTS

Read the result **10 minutes** after dipping. Observations after 11 minutes may be inaccurate due to overdevelopment of the device.

1. Negative result


No line at position T (test): Level of almond, hazelnut, pecan, walnut, cashew and/or pistachio is below the limit of detection. (See **Limitations** section).

Negative

2. Positive result

Any intensity of the red line at position T (test): Level of almond, hazelnut, pecan, walnut, cashew and/or pistachio above detection limit.

Positive

3. Invalid results

If a red line does not appear at position C (control), then the test is invalid.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on below. Training on this product, and all Neogen test kits, is available.

MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts

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