
***Read instructions carefully before starting test***

# **Reveal<sup>®</sup> 3-D**

## **Almond Test**

### **INTRODUCTION AND INTENDED USE**

The Reveal<sup>®</sup> 3-D Almond Test is uniquely designed with 3 lines of detection and can be used virtually anywhere to screen environmental swabs and rinses for the presence of significant levels of almonds. The test's 3-D technology ensures greater reliability with screening than ever before.

The Reveal 3-D Almond Test is intended for use only in an industrial food manufacturing/preparation context or for food labelling enforcement testing. Because of the problems of adequately sampling and extracting almond, it is not suitable for the testing of foods to be consumed in the home or in a restaurant by allergic individuals. The test detects significant (low parts per million, or ppm) levels of almond content in rinses and environmental swabs.

### **DETECTION LIMIT**

The test uses highly specific antibodies to detect almond protein. Almond can be detected from various surfaces using the provided environmental swabs or from rinses.

Utilizing the environmental swabs supplied, levels of 20 µg/ 100 cm<sup>2</sup> of almond on surfaces can be detected.

When analyzing rinses, almond residues are detectable at a level of 5-10 ppm. The presence of cleaners and sanitizers can affect limit of detection in rinses.

### **CROSS REACTIVITY**

The test was used to analyze a panel of potentially cross reacting commodities including nuts, apricot, seeds, legumes and grains. Of these, ground apricot kernel gave a clear cross reactivity with the test. Other samples (Brazil nut, desiccated coconut and pine kernel) gave low positive results but are not considered to be cross contamination (results are possibly due to contamination at source). Further information is detailed in the Reveal 3-D Almond customer validation report, which is available on request.

## **TEST PERFORMANCE**

The Reveal 3-D Almond Test has undergone rigorous validation to evaluate the specificity, sensitivity, robustness and intra- and inter-batch variability of the test method on rinses and environmental swabs. A customer validation report is available on request.

## **SAMPLE COMPATIBILITY**

The Reveal 3-D Almond Test is designed to detect almond on environmental surfaces and in rinses. Although every effort has been made to validate as many variables as possible, there may be some sample types that are not suitable for testing.

Users should perform in-house matrix-specific spike recovery validation work in conjunction with a validated laboratory assay (e.g., Veratox) to help confirm Reveal 3-D Almond Test results. This procedure will highlight any problematic matrices encountered.

The validation of certain food ingredients may be applicable for the Reveal 3-D test format. However, it is recommended that a more broadly inclusive test method, such as Alert or Veratox Almond, be used for food ingredients. Existing commodity validations should be revalidated when suppliers or the manufacturing process has been changed. Please contact a Neogen representative or distributor for additional details.

## **REAGENTS AND MATERIALS**

The Reveal 3-D Almond Test pack contains the following:

1. An instruction leaflet
2. 1 foil pouch, containing 10 yellow Reveal 3-D Almond Test devices
3. 10 sachets containing Type 3 Extraction Buffer
4. 10 sample tubes and caps
5. 10 individually packaged sterile swabs with break off tips

## **SAMPLING TECHNIQUE**

### **1. Clean-in-place rinse water**

As only a small amount of material (0.25 mL) is required for the Reveal 3-D Almond Test, it is important to test a representative portion of liquid.

### **2. Environmental swabbing**

The swabs supplied are intended to be used for the collection of environmental samples from which the presence of almond can be tested. This method can be used to validate the adequacy of cleaning and/or to identify problem areas (e.g., unwanted build-up of almond in processing equipment).

### **3. Food samples**

Contact a Neogen representative or distributor for additional details.

## PRECAUTIONS

1. For environmental testing in an industrial food manufacturing/ preparation, or labeling enforcement context only.
2. Do not use any part of the test beyond the expiry date.
3. Do not open the foil pouch until just before use.
4. Always store the kit between 2–8°C (35–46°F). Avoid freezing.
5. Ensure the foil bag is tightly sealed after removal of a device.

## TEST STORAGE / SAMPLE STABILITY

Store the Reveal 3-D Almond Test kit between 2–8°C (35–46°F). The kit should be used within the expiry date stated on the outer label. Extracted samples should be used within **3 hours** of extraction.

## LIMITATIONS

A negative test from a surface or rinse cannot exclude the possibility that the food contains almond since it may be distributed unevenly on the surface or rinse and may be below the detection limit of the test with that specific sample.


The Reveal 3-D Almond Test is qualitative and should only be used as a preliminary screen for the presence of almond. The validity of results obtained with the test should preferably be viewed in conjunction with data from a quantitative assay.

## RINSE SAMPLING



Remove the following and allow to equilibrate at room temperature before use (**20–30 minutes** out of refrigerator):

- 1 x sample tube
- 1 x Type 3 Extraction Buffer sachet
- 1 x Reveal 3-D Almond Test device (in foil pouch)

1. Carefully tear/cut off the top of the clear Type 3 Extraction Buffer sachet and add the entire contents to the sample tube.
2. Add the 0.25 mL (250 µL) sample to the sample tube. If a pipette is unavailable, 0.25 mL may be crudely estimated by half-filling a white sample tube cap. **NOTE:** Test results will be less accurate.
3. Secure the white cap and shake for **1 minute**.


OR



# SWAB SAMPLING

Remove the following and allow to equilibrate at room temperature before use (20–30 minutes out of refrigerator):

- 1 x sample tube
- 1 x Type 3 Extraction Buffer sachet
- 1 x Reveal 3-D Almond Test device (in foil pouch)
- 1 x sterile swab


1. Carefully tear/cut off the top of the Type 3 Extraction Buffer sachet and add the entire contents to the sample tube.
2. Estimate a swabbing area of approximately 10 cm x 10 cm. Alternatively, use the swab to collect samples of contamination from problem areas (e.g., of processing equipment).
3. Gather the sample with the swab, using one of the following methods:


**For dry surfaces:** Remove a sterile swab from the packaging and wet with extraction buffer in the sample tube. Swab a 10 x 10 cm area by using a crosshatch technique revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing.

**For wet surfaces:** Remove a sterile swab from the packaging and swab a 10 x 10 cm area by using a crosshatch technique revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing. Do not moisten swab prior to use.


4. Return the swab to the Extraction Buffer in the sample tube and carefully break off the moistened end at the prescored mark so it remains in the tube.
5. Secure the cap of the sample tube, taking care to ensure that the stem does not prevent the tube from being properly sealed. Shake for **1 minute**.


## SAMPLE TESTING

1. Remove the lid and fill it with the liquid from the tube. Any froth should remain in the tube.
2. Dip the head of the Reveal 3-D device into the liquid in the lid. Ensure that the cavity is saturated with the liquid.
3. Leave the cavity saturated until you see the liquid running in the test window.
4. Place device on a flat surface and allow test to develop for **5 minutes**.

**NOTE:** It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device.


## READING RESULTS

Liquid will flow into the test window; read the result **5 minutes** after dipping. If distinct lines are not visible at positions O and C, refer to numbers 3 or 4 below.

### 1. Negative result


No line at position T (test): Level of almond below detection limit. (See **Limitations** section).

**Negative**

### 2. Positive result


Any intensity of line at position T (test): Level of almond above detection limit.

**Positive**

### 3. High positive results

No line is visible at position O (overload) and a line is faintly visible or absent at position T: Sample is overloaded with almond.

**High positive**

**High positive**

### 4. Invalid results

If no line appears at position C (control), then the test may be invalid.

## CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

## MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at [www.neogen.com](http://www.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

## WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### TESTING KITS AVAILABLE FROM NEOGEN

#### Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

#### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

#### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

#### Food allergens

- Almonds, crustacea, eggs, gliadin, hazelnut, lupine, milk, mustard, peanut, sesame, soy, walnut

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