



Reveal[®] 3-D Hazelnut Test

For the 'on-site' detection of hazelnut content in food & environmental swabs

The innovative design of the Reveal 3-D Hazelnut gives you the tool to rapidly assess the risks posed by a major food allergen. By identifying cross contamination of ingredient supplies or inadequate cleaning between production batches Reveal 3-D Hazelnut can have a key role as part of a food allergen management system.

Simple & Reliable

The **Reveal 3-D Hazelnut Test** makes it simple to screen for the presence of low levels of Hazelnut in food products and environmental swabs virtually anywhere. Sample preparation and testing take less than **10 minutes**, making it a useful "on-site" screening test. Each kit contains everything needed (including swabs) to perform 10 single use tests.

Increased Sampling & Confidence

Affordability and simplicity of use means the Reveal 3-D test can be implemented more regularly than laboratory food allergen testing. Delivering enhanced confidence and protection through increased sampling and testing of all critical elements of the production process in 'real-time'.

Reveal 3-D Technology

The unique Reveal 3-D technology ensures the three line test will deliver easily recognisable positive / negative results which can be relied upon for real-time decision making, which is crucial when the decision could mean an expensive delay in production or worse, a product recall.

Specificity

No known cross reactivity with a range of nuts, seeds, legumes, and grains. Has strong cross reactivity to walnut (60 ppm whole walnut). On validation the test detected down to 1 ppm in spiked matrices. Please note limit of detection will be matrix specific and validation is recommended for each sample type.

Reveal 3-D Ordering Information



10 tests 902086G 902075M 902082Q 901031P 902087E 901041L 902081S 902093K

Reveal 3-D

WHY?

for real-time decision making

WHEN?

for process control, auditing inspection, troubleshooting and training

WHERE?

in all areas of the process and supply chain

WHO?

by anyone, anywhere, no training or equipment required

BEST PRACTICE

Reveal 3-D tests are uniquely designed for on-site validation and verification of allergen controls because best practice requires specific detection of food allergen proteins

Make Reveal 3-D an integral part of your Allergen Control Plan

FOR FURTHER INFORMATION:

Neogen Corporation

US: 800/234-5333 • 517/372-9200
Email: foodsafety@neogen.com
Website: www.neogen.com

Neogen Europe Corporation

UK: +44 1292 525 275
Email: customerservice@neogeneurope.com
Website: www.neogeneurope.com

